


# AviAir

Single Use Bougies



# AviAir

## Single Use Bougies



Tracheal intubation is difficult in 1–4% and impossible in 0.05–0.35% of patients who have seemingly normal airways. <sup>[1]</sup>

The use of a tracheal tube introducer or bougie early on in such cases is recommended. <sup>[2]</sup>

AviAir is a range of single use bougies which have been developed with the benefits of traditional gum elastic bougies in mind.

AviAir bougies have been designed with clear graduated markings and a luminescent coude tip to assist with placement and visualisation during clinical use.

### Key Features\*:

- Single Use
- Luminescent coude tip\*\*
- Clear graduated markings on left hand side
- Superior memory
- Good flexibility



[1] Benumof JL. Management of the difficult adult airway. With special emphasis on awake tracheal intubation. *Anesthesiology* 1991; 75: 1087–110.

[2] Henderson JJ, Popat MT, Latto IP, Pearce AC. Difficult Airway Society guidelines for management of the unanticipated difficult intubation. *Anaesthesia* 2004; 59: 675–94.

### Product Information

Code	Description	Box Qty
AMTI10800	Intubating bougie 10ch length 800mm with soft luminescent tip - Coming Soon	10
AMTI15750	Intubating bougie 15ch length 750mm with soft luminescent tip	10
AMTI15750V	Intubating bougie, vented 15ch length 750mm - Coming Soon	10
INBV147500X	Introducer with angled (coude) tip, vented with O <sub>2</sub> connector, 14ch 750mm	10

\* Key Features relate to AMTI Product Codes

\*\* Non-vented versions



Distributed by:  
Armstrong Medical  
Wattstown Business Park  
Newbridge Road  
Coleraine BT52 1BS  
Northern Ireland

T +44 (0) 28 7035 6029

F +44 (0) 28 7035 6875

E [info@armstrongmedical.net](mailto:info@armstrongmedical.net)

W [www.armstrongmedical.net](http://www.armstrongmedical.net)

Creating Support for Life